

HealthLab Candy-Cable VSS-03



Technical Data

Designation

HealthLab Bus Connection Cable
Master / Satellite / Satellite

Cable Diameters

2 x AWG-26 power supply
2 x AWG-28 data line

Climatic Conditions

acc. to DIN EN 60204-1 (10-2014)

Ambient Temperature

operation: -20 ... +55 °C
transport / storage: -25 ... +60 °C

Humidity / Altitude

20 ... 90 % RH (without condensing)
up to 8.000 m

Dimensions

wire length: approx. 300 mm
weight: approx. 8 g

Connection Cable HealthLab Serial Slave Bus

By means of the HealthLab Candy-Cable VSS-03, the system bus connection (HealthLab Serial Slave Bus) between HealthLab Master (e.g. HL Satellite-Master SAT-0M) and HealthLab Measuring-Satellite, as well as between Satellite and Satellite is done. Both, the power supply of the connected devices, as well as the data exchange between master and satellites takes place via this system bus connection cable.



Note: The psycho-physiological monitoring system HealthLab is manufactured and delivered in configurations according to customer's request. The HealthLab components are not certified for use in the medical field. Therefore they shall be used solely for research purposes in scientific area.

Links on further documents:

- Hardware: Master HFM-01 → <https://secure.turboj.de/documents/HFM-01.pdf>
- Software: Heally Control → https://secure.turboj.de/documents/Heally5_en.pdf

Ordering Information:

VSS-03.Candy-Cable:
(system bus connection
cable Master / SAT / SAT)

Part-No.

E1816